



A Symbol of Quality



CENTRIFUGE



AIISHIL INTERNATIONAL

18/19/20, Shreeji Arcade, 2nd Floor, Opp. Nitin Co. Thane (w) 400602. INDIA.

Phone : 022-25390474, 25390475, Fax : 022-25390476

E-mail : infob@aiishilgroup.com

CENTRIFUGE

"AIISHIL" CENTRIFUGE MACHINES (MEDICO/CLINICAL):

JEQ-6 (a): Dome type, sturdy, stable, heavy body with silently mounted motor for vibration-free performance with builtin 5-speed regulator and one non-removeable swingout rotor. Complete with tube carriers. Maximum speed 3500 RPM. Neatly finished in stoving enamel paint.



JEQ-6 (a)

Medical/Clinical	(i)	(ii)	(HI)
Capacity (ml.)	4x15	6 x 15	8x15
Max. RCF'g'	1850	1600	1600
Doctor			
Max. RCF'g"	1250	1250	1250



JEQ-6(C)

"AIISHIL" COMPACT LABORATORY CENTRIFUGE :

JEQ-6(C) :Compact rectangular body made of mild steel duly powder coated, with clear acrylic lockable lid, continuously variable speed regulator, safety fuse, pilot lamp etc. provided on streamlined panel. With plug-cord to work on 220 V AC 50 Hz single phase. Maximum speed 4000 RPM. A special safety mode is provided that when lid is open instrument will stop working.

Capacity (ml.) 4 x 50 6 x 50 8 x 15 12 x 15 16 x 15

Optionals :

1. Automatic 0-60 minutes Timer.
2. Speedometer Digital.
3. Maximum speed 10000 can also be provided.

CENTRIFUGE

AIISHIL" REFRIGERATED CENTRIFUGE :

JEQ-6 (D) : Streamlined rectangular body with elegant panel and clear acrylic lockable lid. Complete with pilot lamp, continuously variable speed regulator, Digital speed-indicator, plug-cord to work on 220 V AC 50 Hz single phase. A hermetically sealed compressor for cooling from ambient to -200C. Temperature controlled by Digital Temp. Controller cum Indicator with accuracy $\pm 0.50C$. Maximum speed 5000 RPM. The centrifuge has an electrical interlock to ensure that the unit can not run with cover open.



JEQ-6 (D)

- (a) Capacity : 8x15 ml. (b) Capacity : 12x15 ml.
(c) Capacity : 16x15 ml. (d) Capacity : 12x1.5 ml.
(e) Capacity : 16x1.5 ml. (f) Capacity : 8x50 ml.

Note : Use of suitable Voltage stablizer is a must.

Optional : Automatic Voltage Stabilizer 2 KVA.

OUR OTHER PRODUCTS RANGE



☞ VISUAL EQUIPMENT (PROJECTORS).



☞ MICROTOMY EQUIPMENTS.

☞ MICROSCOPY EQUIPMENTS.

☞ CLEAN AIR SYSTEM.

☞ HEAT & REFRIGERATION SYSTEM.

☞ BLOOD BANKING.



☞ GENERAL LABORATORY SYSTEM.

☞ PHARMACEUTICAL & PHARMACOLOGY EQUIPMENTS.



☞ MEDICAL ANATOMY INSTRUMENTS.

☞ ELECTRONIC ANALYTICAL INSTRUMENTS.

☞ INCINERATOR.



☞ ENGINEERING COLLEGE EQUIPMENTS.



AIISHIL INTERNATIONAL

18/19/20, Shreeji Arcade, 2nd Floor, Opp. Nitin Co. Thane (w) 400602. INDIA.

Phone : 022- 25390474, 5390475, Fax : 022-25390476

E-mail : infob@aishilgroup.com