



A Symbol of Quality

# MEDICAL ANATOMY INSTRUMENTS



BONE & MEAT CUTTING MACHINE



MORTUARY CHAMBER



AUTOPSY TABLE



COLD ROOM



CADAVEROUS INJECTOR



## AIISHIL INTERNATIONAL

18/19/20, Shreeji Arcade, 2nd Floor, Opp. Nitin Co. Thane (w) 400602, INDIA  
Phone : 022-25390474, 25390475, Fax : 022-25390476  
E-mail : [infob@aiishilgroup.com](mailto:infob@aiishilgroup.com)

## "AIISHIL" MORTUARY CHAMBER :

Inner Chamber and tray size 200x56x8 cms. sliding on stainless steel U-rails are made of stainless steel. Outer chamber made of mild steel duly painted. Temperature range - 5°C to ambient controlled through digital temperature controller-cum-indicator. Complete with built in stabilizer & defrosting system. Chamber is provided with locked door & has index card holder & magnetic gasket. Workable on 220V AC 50 Hz single phase. Model for single body, two bodies, three bodies, four bodies & six bodies are available.

**Models available : JKA-081, 082, 083, 084, & 086**



**JKA - 081**

## "AIISHIL" BONE & MEAT CUTTING MACHINE :

**JEQ-11** : Used by Anatomy & Meat departments for cutting & packing of meat, fish etc. Sturdy base with large moving and extensive table working on ball-bearing rollers. Table made of aluminium or stainless steel sheet.

Supplied with one blade, starter, cord-plug etc. to work on 220 VAC 50 Hz single phase.

Table size : 785x585 mm

Table travel : 1245 mm.

**Optional** : Spare cutting blade.



**JEQ-11**

## "AIISHIL" MEAT CUTTING MACHINE / CAUTRY MACHINE :

**JEQ : 11B** : Constructed on heavy cast iron angles covered with thick gauge mild steel sheet duly powder coated.

Electrically operated.

Stainless Steel Top.

Adjustment of the cutter for upward & downward position.

Total unit on casters for easy mobility.

## "AIISHIL" AUTOPSY TABLE

**JEQ : 11T** : Autopsy stainless steel table on S.S. tubular stand of first class workshop scientifically designed for maximum efficiency. Two nos. 6" wide body, support, Head rest, Crusher for disposal swing spout mixing faucet. A Stainless Steel Water basin with Hot & Cold mixing tap. Complete with Drain Pipe.



**JEQ-11T**

## "AIISHIL" COLD ROOM

**JEQ : 11CR** : Total Body with Heavy sturdy Angle Iron & Channels outer body is made of MS sheet duly painted. Inner chamber made S.S. Single door opening system with Door Opening Alarm. Temperature 4°C is brought down by Hermetically sealed compressor and controlled by microprocessor based temperature variation alarm facility controller, Puf Insulation to maintain the temperature, inner chamber illumination total unit is moveable on Casters. Capacity 108 Sq. Feet.



**JEQ-11CR**

## "AIISHIL" METAL MINCING MACHINE :

**JEQ-31** : The unit is mounted on Angle iron to give the sturdiness to the instruments it is ment for crushing and destroying the used small blades, needles and accessories.

## "AIISHIL" CADAVEROUS INJECTOR

**JEQ-41** : For injecting Formaldehyde solution in Cadaverous at much higher speed than normal gravity process. Unit is fully covered & mounted on a portable trolley having four castor wheels for easy movement.

Unit consists of one air compressor fitted with 1/2 hp Motor which is connected with a Stainless Steel Tank of 10 liters capacity meant for storing and injecting the solution.

Tank is fitted with a safety valve, pressure gauge and rubber tubing having provision for injection.

Supplied complete with electric cord, plug and suitable to work on 220 volts, 1ph 50 Hz, Ac supply.



**JEQ-41**

# OUR OTHER PRODUCT RANGE :



VISUAL EQUIPMENTS

MICROTOMY EQUIPMENTS

MICROSCOPY EQUIPMENT



CLEAN AIR SYSTEM

HEAT & REFRIGERATION  
SYSTEM



BLOOD BANKING

GENERAL LABORATORY  
SYSTEM



PHARMACEUTICAL

PHARMACOLOGY EQUIPMENTS



ELECTRONIC ANALYTICAL  
INSTRUMENTS

INCINERATORS

ENGINEERING COLLEGE  
EQUIPMENTS



## AIISHIL INTERNATIONAL

18/19/20, Shreeji Arcade, 2nd Floor, Opp. Nitin Co. Thane (w) 400602. INDIA.  
Phone : 022-25390474, 25390475, Fax : 022-25390476  
E-mail : [infob@aishilgroup.com](mailto:infob@aishilgroup.com)